- EN** more information on [www.fkg.ch](http://www.fkg.ch)
- FR** informations supplémentaires sur [www.fkg.ch](http://www.fkg.ch)
- DE** weitere Informationen auf [www.fkg.ch](http://www.fkg.ch)


FKG REF. - 99.7AA.10.14A.XX - notice n°103 - 2013/09


## BioRace

Description - Beschreibung

### BioRace Basic Set

### BioRace Extended Set

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## EN Golden rules

Speed: 600 rpm - Torque: 1 Ncm

- > Max. 3-4 back and forth gentle strokes
- > Light touch, let the instrument work
- > Work 3-4 seconds in a row, get out, clean the blade and irrigate the canal

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### Prepare glide path first, then

Fill the canals and pulp chamber with irrigant.

1. BR0 - only 4 gentle strokes - clean the flutes.
  - Repeat until approximately 4-6mm of coronal part of the canal has been prepared.
  - After use of BR0, repeat irrigation.
  - Recapitulate to full working length (WL) with a manual file #15.
  - Fill the canal and pulp chamber with irrigant.
2. Use BR1 with 4 gentle strokes. If this instrument does not reach the WL, clean the instrument and repeat until the WL is achieved.
3. Use BR2 and BR3 as described for BR1. DO NOT use BR3 to full WL on canals with severe apical curvatures. In this case go to step 1 of BioRace Extended Set protocol.
4. Use BR4 and BR5 as explained for BR1-3. In most cases, the final apical preparation is then achieved.
5. For larger canals, (see anatomical chart on [www.fkg.ch](http://www.fkg.ch)), go to step 3 of BioRace Extended Set protocol.

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## BioRace Extended Set

1 blister includes/contient/enthält: 1 × BR4C, 1 × BR5C, 1 × BR6, 1 × BR7

1. For severe apical curvatures, instruments BR4C and BR5C should be used to prepare the apical canal. If the instrument does not reach the WL with 4 gentle strokes, DO NOT FORCE the instrument. Irrigate the canals and repeat.
2. For complicated curvatures it is recommended to use additional FKG instruments (e.g. S-Apex, ScoutRace or manual SMG files).
3. For larger canals the two additional instruments BR6 and BR7 from BioRace Extended Set can be used as explained for BR1-5.

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NB. Copious irrigation at all times. Cleaning the files after 4 gentle strokes is essential for safe and efficient use of these instruments.

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EN

## Sterilization protocol for reusable medical devices

### Non sterile products

#### First use

Carefully follow the rules herebelow  
**Cleaning and Sterilization**

### Sterile products

#### First use

1. Open the alveolus
2. Extract the instrument (use of gloves imperative)
3. Use the instrument according to good dental practice

#### Second use:

Carefully follow the rules herebelow  
**Second and following uses.**

### Cleaning

1. Manual cleaning with/without Ultrasonic assistance  
Thorough cleaning allows efficient disinfection & sterilization.
2. Rinse  
Under distilled/demineralized water at least 1 min.  
Then dry the instruments.

3. Inspection - Check  
Sort out damaged or worn out instruments.

4. Packing  
Place the instruments in a suitable device/container and pack the devices in sterilisation pouches ISO 11607-1.

### Sterilization

5. According EN/ISO 17664 protocol  
Autoclave: 134°C/273°F,  
2.2 bar during 18 min.

6. Storage  
Keep devices in sterilization wrap/pouch in a dry and clean environment.

 Always refer to manufacturer's instructions for cycle & duration

### Second and following uses

7. Pre-disinfection  
Plunge instruments after use in detergent solution and eventually brush them manually

8. Rinse  
Under running water at least 1 min.  
Then dry the instruments.
9. Follow steps 1 - 6 above

 Complete operating protocol and warnings on [www.fkg.ch](http://www.fkg.ch)

### General information

#### Pre-/Disinfection DO NOT USE

- Solutions containing:
- phenol (corrosion)
  - aldehyde (blood fixation)
  - di-/triethanolamines (corrosion)

**Sterilization with chemiclav or dry air devices** have not been validated by FKG Dentaire.  
Refer to manufacturer's instructions for cycle & duration. Nevertheless our instruments withstand such methods.

#### Caution:

This product contains nickel and should not be used for individuals with known allergic sensitivity to this metal.

A single patient use is recommended to avoid cross-contamination.

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