



Instruments for Injection




ASPIJECT® Active

Injection syringe for infiltration
and block anaesthesia

ASPIJECT® Active is an active aspirating syringe using a standard dental cartridge. The piston end is supplied with a self-tapping spiral-hook which engages with the plunger. Manufactured in stainless steel to the highest of engineering standards. Its unique styling and design offers the ultimate in balance, feel and tactile control.

- Made from high-grade stainless steel and heat-resistant engineering polymer
- Supplied with a self-tapping spiral-hook which engages with the plunger
- New superior > no friction < surface treatment of the piston
- Improved design with contra handle for optimal control
- Can be sterilized by any method
- Material and construction guaranteed for 5 years

ASPIJECT® Active



RØNVIG
Dental Mfg. A/S

Gl. Vejlevej 59 · DK-8721 Daugaard · Denmark
Tel.: (+45) 70 23 34 11 · Fax: (+45) 76 74 07 98
e-mail: export@ronvig.com · www.ronvig.com