



# MixStar-eMotion

**High-Speed Design.**

HIGH-SPEED DESIGN.



2012  
**REALITY**  
Four Star Award  
★★★★

## Fully automatic mixing with electronic perfection.

Mixing has never been easier, faster and so safe: With the MixStar-eMotion DMG presents the newest generation of fully automatic impression material mixing devices.

This high-tech unit is extremely easy to operate. The MixStar-eMotion can be controlled with only one hand by means of a quick push-button operation – incorrect dosages are impossible. And its stylish design impresses optically as well.

## Versatility included.

Regardless of your impression material and technique preferences: the MixStar-eMotion provides true relief for all daily practice routines. For fully automatic applications DMG offers a comprehensive line of high-quality impression materials. And the device is compatible with commercially available, pre-filled cartridges and foil bags.

| Features  | Benefits   |
|---|--|
| Fully automatic   | Extremely easy to operate<br>Hand-free operation                                 |
| 3 speed levels  | Extremely fast and time efficient<br>Tailored to the dental team's varying needs |
| 2 timers for working and setting time                   | Error-free treatment process   |
| Invariably homogeneous and bubble-free material quality | Optimal dosage<br>Perfect mixing ratio   |
| USB interface to update the software                    | Always current   |



Honigum and StatusBlue MixStar cartridges for dynamic mixing in the MixStar-eMotion.



Easy loading of a syringe at speed level 1.



Quick loading of a tray at speed level 2 or 3.

# HIGH QUALITY FOR ALL APPLICATIONS.

## Honigum: An entire family for highest precision

Typically, impression materials are either especially flowable or stable. Honigum, however, overcomes these opposites. Thanks to its patented rheologically active matrix Honigum yields best results in both disciplines. For this and other reasons Honigum received the best clinical rating\* of all 50 VPS tested by the noted »DENTAL ADVISOR«. In addition, due to its excellent fixation characteristics Honigum-Heavy is recommended for implant impressions by the prestigious Brånemark Center. And with Honigum-MixStar Putty DMG offers a dynamically mixed true putty material in the MixStar cartridge.

## StatusBlue: The VPS for all alginate indications

Due to its high dimensional stability, impressions taken with StatusBlue neither shrink nor swell thus allowing for easy long-term storage of the impression for weeks and ensuring exact fit.

The material's storage stability also makes multiple casts of the impression possible. Repeat or multiple impressions, e.g. due to a broken provisional, are therefore unnecessary. In a comparison on the topic »Alginates and Alginate Substitutes« conducted by the US test institute »THE DENTAL ADVISOR« StatusBlue scored the highest ratings of all materials.\*\*

Besides StatusBlue and Honigum, the tried and tested Silagum-MixStar Putty Soft is available for use with the MixStar-eMotion as well.

**The optimized MixStar-cartridges** are color-coded for a quick product identification. In addition, the labels indicate the respective working and setting times for easier orientation and in order to prevent errors.

\*The Dental Advisor, Vol. 23, No. 3

\*\*The Dental Advisor, Vol. 24, No. 2, p. 2-7



Color-coded locking-systems and labels indicate the respective working times.



Detailed reproduction of the oral situation with Honigum-MixStar Putty and Honigum-Light.



Transfer impression with Honigum-Heavy.



Exact situation impression with StatusBlue.



# STATE-OF-THE-ART AUTOMATIC MIXING.

## Easier.

The MixStar-eMotion is extremely easy to operate. One quick push-button operation and the material is automatically and precisely dosed, mixed and dispensed. And, compared to other devices available on the market, the Start button does not have to be held during the mixing and delivery process thus allowing for a very comfortable and hands-free workflow. For additional relief the unit features two timer functions for the working and setting times which allow you to program, store and easily retrieve the respective times for the different impression materials.

## Faster.

The MixStar-eMotion is the only mixing device featuring three different speed levels, the highest of which delivering twice as fast as its predecessor model, the MixStar. For loading the tray, for example, you can save up to 20 seconds thus gaining additional working time for the material. The different speed levels were implemented in order to optimize the impression taking process. Whether you load syringes or trays, you're an experienced or novice assistant: The MixStar-eMotion offers the right speed for every situation.

## Safe.

In an electronically controlled manner the MixStar-eMotion provides an invariably homogeneous and bubble-free material quality. The material always comes optimally dosed into the tray or syringe – and you can be sure to have worked economically since material waste is reduced to a minimum. By storing the correct working and setting times for the individual materials mixing and impression errors can be prevented.

The DMG MixStar cartridges indicate the respective times in order to support an error-free treatment process.



Display with 2 timers and 3 different speed levels.



Button for a quick cartridge exchange.



Start and Stop button for automatic mixing.



Connector for the optional foot pedal and USB interface to update the software.

## PACKAGING

### MixStar-eMotion

#### Package

|                   |            |
|-------------------|------------|
| 1 MixStar-eMotion |            |
| 230 V int.        | REF 900480 |
| 230 V GB          | REF 900483 |
| 230 V NZ/AUS      | REF 900484 |
| 115 V int.        | REF 900482 |

#### Accessories

|                 |            |
|-----------------|------------|
| 50 MixStar-Tips | REF 909567 |
| 1 Foot pedal    | REF 900310 |

### Honigum-Light

#### Package

|                                    |            |
|------------------------------------|------------|
| 4 Cartridges @ 25 ml               |            |
| 16 Automix-Tips, 16 Intraoral-Tips |            |
| Honigum-Light                      | REF 909539 |
| Honigum-Light Fast                 | REF 909639 |

|                                    |            |
|------------------------------------|------------|
| 2 Cartridges @ 50 ml               |            |
| 12 Automix-Tips, 12 Intraoral-Tips |            |
| Honigum-Light                      | REF 909831 |
| Honigum-Light Fast                 | REF 909832 |

#### Economy Pack

|                                    |            |
|------------------------------------|------------|
| 16 Cartridges @ 25 ml              |            |
| 50 Automix-Tips, 50 Intraoral-Tips |            |
| Honigum-Light                      | REF 909690 |
| Honigum-Light Fast                 | REF 909691 |

|                                    |            |
|------------------------------------|------------|
| 8 Cartridges @ 50 ml               |            |
| 50 Automix-Tips, 50 Intraoral-Tips |            |
| Honigum-Light                      | REF 909835 |
| Honigum-Light Fast                 | REF 909836 |

### Honigum-MixStar Putty

#### Package

|                       |            |
|-----------------------|------------|
| 1 Cartridges @ 380 ml |            |
| 10 MixStar-Tips       | REF 919374 |

#### Economy Pack

|                       |            |
|-----------------------|------------|
| 5 Cartridges @ 380 ml |            |
| 50 MixStar-Tips       | REF 919375 |

### Honigum-MixStar Heavy

#### Package

|                                       |            |
|---------------------------------------|------------|
| 1 Cartridge @ 380 ml, 10 MixStar-Tips |            |
| Honigum-Heavy                         | REF 909537 |
| Honigum-Heavy Fast                    | REF 909638 |

#### Economy Pack

|  |            |
|--|------------|
| 5 Cartridges @ 380 ml, 50 MixStar-Tips |            |
| Honigum-Heavy                          | REF 909538 |
| Honigum-Heavy Fast                     | REF 909637 |

### Honigum-MixStar Mono

#### Package

|                      |            |
|----------------------|------------|
| 1 Cartridge @ 380 ml |            |
| 10 MixStar-Tips      | REF 909568 |

#### Economy Pack

|                       |            |
|-----------------------|------------|
| 5 Cartridges @ 380 ml |            |
| 50 MixStar-Tips       | REF 909569 |

### StatusBlue

#### Package

|                      |            |
|----------------------|------------|
| 1 Cartridge @ 380 ml |            |
| 10 MixStar-Tips      | REF 909693 |

#### Economy Pack

|                       |            |
|-----------------------|------------|
| 5 Cartridges @ 380 ml |            |
| 50 MixStar-Tips       | REF 909694 |



## DMG

### Chemisch-Pharmazeutische Fabrik GmbH

Elbgaustraße 248 22547 Hamburg Germany

Fon: +49. (0) 40. 84 006-0

Fax: +49. (0) 40. 84 006-222

info@dmg-dental.com

www.dmg-dental.com